

ANDOCOR



CARDIOVASCULAR CANNULATION PRODUCTS	3 - 9
Arterial cannulae	3
Aortic catheter	7
Flex line venous cannulae	8
Flex line two stage venous cannulae	9
CARDIOPLEGIA PRODUCTS	10 - 16
Aortic root cannulae	10
Retrograde cardioplegia cannula	12
Cardioplegia needle	12
Cardioplegia adapters	13
Ostial perfusion cannulae	16
VENT CATHETERS	17 - 20
Silicone vent catheter	17
PVC vent catheters	18
SUCKERS	21 - 23
Rigid suckers	21
Flexible suckers	21
Yankauer suction tube	22
Pericardial sump	23
GAS DIFFUSER	24
HEMOCONCENTRATORS	25
ACCESSORIES	26 - 29
Suction connecting tubes	26
Aspiration tubes	26
Porous vent plugs	26
Vented connector cap	26
Vessel cannulae	27
Tubing organizer	27
Connectors	28
Cannulation tourniquet sets	29
Quick prime line	29
Pressure monitoring line	29
MAMMARIA RETRACTOR	30

CARDIOVASCULAR CANNULATION PRODUCTS



REINFORCED, CURVED TIP WITH FLANGE

ARTERIAL CANNULAE

AndoCor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant reinforced tubing
- ✓ Easy tip orientation thanks to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP
WITHOUT CONNECTOR	20 Fr	33cm	A2016	A2016V		24 cm	B2016	B2016V	
	22 Fr	33cm	A2216	A2216V		24 cm	B2216	B2216V	
	24 Fr	33cm	A2416	A2416V		24 cm	B2416	B2416V	
WITH CONNECTOR	20 Fr	33cm	A20161		A20161V	24 cm	B20161		B20161V
	22 Fr	33cm	A22161		A22161V	24 cm	B22161		B22161V
	24 Fr	33cm	A24161		A24161V	24 cm	B24161		B24161V
WITH LUER CONNECTOR	20 Fr	33cm	A20162		A20162V	24 cm	B20162		B20162V
	22 Fr	33cm	A22162		A22162V	24 cm	B22162		B22162V
	24 Fr	33cm	A24162		A24162V	24 cm	B24162		B24162V

CARDIOVASCULAR CANNULATION PRODUCTS



REINFORCED, STRAIGHT TIP WITH RING

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant reinforced tubing
- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP	LENGTH	WITHOUT VENT PLUG/ RED CAP	WITH VENT PLUG	WITH VENTED RED CAP
WITHOUT CONNECTOR	20 Fr	33 cm	A2013	A2013V		24 cm	B2013	B2013V	
	22 Fr	33 cm	A2213	A2213V		24 cm	B2213	B2213V	
	24 Fr	33 cm	A2413	A2413V		24 cm	B2413	B2413V	
WITH CONNECTOR	20 Fr	33 cm	A20131		A20131V	24 cm	B20131		B20131V
	22 Fr	33 cm	A22131		A22131V	24 cm	B22131		B22131V
	24 Fr	33 cm	A24131		A24131V	24 cm	B24131		B24131V
WITH LUER CONNECTOR	20 Fr	33 cm	A20132		A20132V	24 cm	B20132		B20132V
	22 Fr	33 cm	A22132		A22132V	24 cm	B22132		B22132V
	24 Fr	33 cm	A24132		A24132V	24 cm	B24132		B24132V

CARDIOVASCULAR CANNULATION PRODUCTS



NON-REINFORCED, CURVED TIP WITH FLANGE

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant pvc tubing
- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENTED RED CAP	WITH VENTED RED CAP
WITHOUT CONNECTOR	20 Fr	25 cm	A2006	
	22 Fr	25 cm	A2206	
	24 Fr	25 cm	A2406	
WITH CONNECTOR	20 Fr	25 cm	A20061	A20061V
	22 Fr	25 cm	A22061	A22061V
	24 Fr	25 cm	A24061	A24061V
WITH LUER CONNECTOR	20 Fr	25 cm	A20062	A20062V
	22 Fr	25 cm	A22062	A22062V
	24 Fr	25 cm	A24062	A24062V

CARDIOVASCULAR CANNULATION PRODUCTS



NON-REINFORCED, STRAIGHT TIP WITH RING

ARTERIAL CANNULAE

Andocor arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS	
Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant pvc tubing
- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth plastic tip



	SIZE	LENGTH	WITHOUT VENTED RED CAP	WITH VENTED RED CAP
WITHOUT CONNECTOR	20 Fr	25 cm	A2003	
	22 Fr	25 cm	A2203	
	24 Fr	25 cm	A2403	
WITH CONNECTOR	20 Fr	25 cm	A20031	A20031V
	22 Fr	25 cm	A22031	A22031V
	24 Fr	25 cm	A24031	A24031V
WITH LUER CONNECTOR	20 Fr	25 cm	A20032	A20032V
	22 Fr	25 cm	A22032	A22032V
	24 Fr	25 cm	A24032	A24032V



AORTIC CATHETER

Andocor aortic catheters are intended for perfusion of the ascending aorta during cardiopulmonary bypass.




DETAILS

Connection site	3/8"
Qty/box	20 Pcs



Features & Benefits

- ✓ Easy tip orientation due to indexing marker
- ✓ Smooth beveled tip
- ✓ Movable suture ring




	SIZE	LENGTH	WITHOUT VENTED RED CAP	WITH VENTED RED CAP
WITH CONNECTOR	20 Fr	24 cm	AC201341	AC201341V
	22 Fr	24 cm	AC221341	AC221341V
	24 Fr	24 cm	AC241341	AC241341V
	26 Fr	24 cm	AC261341	AC261341V
	30 Fr	24 cm	AC301341	AC301341V
WITH LUER CONNECTOR	22 Fr	24 cm	AC221342	AC221342V
	24 Fr	24 cm	AC241342	AC241342V
	26 Fr	24 cm	AC261342	AC261342V
	30 Fr	24 cm	AC301342	AC301342V

SPECIAL CONFIGURATION: HARD TIP

WITH CONNECTOR	24 Fr	24 cm	AC241341H	AC241341VH
WITH LUER CONNECTOR	24 Fr	24 cm	AC241342H	AC241342VH

CARDIOVASCULAR CANNULATION PRODUCTS



REINFORCED WITH LIGHTHOUSE TIP


FLEX LINE VENOUS CANNULAE

Andocor venous cannulae are intended to drain venous blood from the right atrium and caval veins.

DETAILS	
Length	39 cm
Qty/box	20 Pcs

+ Features & Benefits

- ✓ Single piece construction
- ✓ Open and flexible lighthouse tip
- ✓ Easy positioning
- ✓ Insertion depth marker



SIZE	STRAIGHT TIP	RIGHT ANGLED 90° TIP	CONNECTION SITE
20 Fr	01V201L8	01V201L9	3/8"
22 Fr	01V221L8	01V221L9	3/8"
24 Fr	01V241L8	01V241L9	3/8"
26 Fr	01V261L8	01V261L9	3/8"
28 Fr	01V281L8	01V281L9	3/8"
30 Fr	01V301L8	01V301L9	3/8"
32 Fr	01V321L8	01V321L9	3/8"
34 Fr	01V341L8	01V341L9	3/8"
36 Fr	01V361L8	01V361L9	1/2"

SIZE	NON-REINFORCED	CONNECTION SITE
32 Fr	VC3208	3/8"



PROXIMAL AND DISTAL REINFORCED WITH LIGHTHOUSE TIP

FLEX LINE TWO STAGE VENOUS CANNULAE


Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.




Features & Benefits

- ✓ Single piece construction
- ✓ Open and flexible lighthouse tip
- ✓ Easy positioning
- ✓ Obturator
- ✓ Insertion depth marker

DETAILS	
Connection site	1/2"
Length	40 cm
Qty/box	20 Pcs



	SIZE	STANDARD BODY	FLAT BODY
WITHOUT CONNECTOR	32/40 Fr	01V32L40	01V32L40FB
	34/46 Fr	01V34L46	01V34L46FB
	36/46 Fr	01V36L46	01V36L46FB
	36/51 Fr	01V36L51	01V36L51FB
WITH CONNECTOR	32/40 Fr	01V32L401	01V32L401FB
	34/46 Fr	01V34L461	01V34L461FB
	36/46 Fr	01V36L461	01V36L461FB
	36/51 Fr	01V36L511	01V36L511FB
WITH LUER CONNECTOR	32/40 Fr	01V32L402	01V32L402FB
	34/46 Fr	01V34L462	01V34L462FB
	36/46 Fr	01V36L462	01V36L462FB
	36/51 Fr	01V36L512	01V36L512FB




DETAILS	
Length	15 cm
Qty/box	20 Pcs

WITH OR WITHOUT VENT LINE

AORTIC ROOT CANNULAE

Andocor aortic root cannulae are intended for use during cardiopulmonary bypass for antegrade delivery of cardioplegia solutions and venting of the heart.

The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.



	CODE	TIP SIZE	NEEDLE SIZE
WITHOUT VENT LINE	ARC07	6 Fr	16 G
	ARC04	7 Fr	14 G
	ARC01	9 Fr	12 G
WITH VENT LINE	ARC08	6 Fr	16 G
	ARC05	7 Fr	14 G
	ARC02	9 Fr	12 G

CARDIOPLEGIA PRODUCTS



WITH ANTI-BLEEDING VALVE



Features & Benefits

✓ Ease of use



	CODE	TIP SIZE	NEEDLE SIZE
WITHOUT VENT LINE	ARC07V	6 Fr	16 G
	ARC04V	7 Fr	14 G
	ARC01V	9 Fr	12 G
WITH VENT LINE	ARC08V	6 Fr	16 G
	ARC05V	7 Fr	14 G
	ARC02V	9 Fr	12 G

SPECIAL CONFIGURATIONS




	CODE	TIP SIZE	NEEDLE SIZE	LENGTH
LONG TYPE WITHOUT VENT LINE	ARC14	7 Fr	14 G	28 cm
WITH VENT LINE WITH 2 CLAMPS	ARC10	7 Fr	14 G	15 cm



RETROGRADE CARDIOPLEGIA CANNULA

Andocor retrograde cardioplegia cannulae are intended to be connected to the cardioplegia line in order to infuse cardioplegic solution and blood into the patient's coronary sinus during open heart procedures.



Features & Benefits

- ✓ Self-inflating
- ✓ Silicone balloon
- ✓ Malleable stylet

DETAILS	
Code	RCC01
Length	31 cm
Balloon diameter	18 mm
Size	14 Fr
Qty/box	20 Pcs

CARDIOPLEGIA NEEDLE

Andocor cardioplegia needles are intended for use in delivery of cardioplegia solution during cardiopulmonary bypass surgery.



Features & Benefits

- ✓ Protective cap
- ✓ Stainless steel needle
- ✓ Luer lock connection

DETAILS	
Code	CN16
Length needle	1,4 cm
Size needle	16G
Total length	5 cm
Qty/box	20 Pcs

CARDIOPLEGIA ADAPTERS

Andocor cardioplegia adapters connect to the aortic root cannulae or vessel cannulae and are intended for the delivery of cardioplegia solution or venting the heart during cardiopulmonary bypass.

DETAILS	
Qty/box	20 Pcs

DIFFERENT CONFIGURATIONS

CSY03

PERFUSION ADAPTER WITH 2 MALE LUER LOCK AND 1 FEMALE LUER LOCK CONNECTOR AND COLOR-CODED CLAMPS



CSY04

MULTIPLE PERFUSION ADAPTER FOR AORTIC ROOT AND THREE VEIN GRAFTS



CSY05

MULTIPLE PERFUSION ADAPTER FOR AORTIC ROOT AND THREE VEIN GRAFTS WITH ATTACHED VESSEL CANNULAE



CSY06

DUAL PERFUSION SET WITH MALE LUER LOCK, CONE CONNECTOR AND COLOR CODED CLAMPS



CSY07

MULTIPLE PERFUSION ADAPTER WITH VENT LINE AND ATTACHED VESSEL CANNULAE



CSY08

AORTIC PERFUSION SET WITH VENT LINE, TOTAL LENGTH = 30 CM



CSY09

DUAL PERFUSION SET WITH MALE LUER LOCK, FEMALE LUER LOCK CONNECTORS AND COLOR-CODED CLAMPS



CSY10

MULTIPLE SET FOR PERFUSION OF THREE VEIN GRAFTS



CSY11

MULTIPLE SET FOR PERFUSION OF THREE VEIN GRAFTS WITH ATTACHED VESSEL CANNULAE



CSY12

MULTIPLE PERFUSION ADAPTER FOR RETROGRADE AND AORTIC ROOT WITH VENT LINE



CSY13

AORTIC PERFUSION SET WITH VENT LINE, TOTAL LENGTH = 18 CM



CSY14

AORTIC PERFUSION SET WITH 1/4"x1/4"x1/4" CONNECTOR



CSY15

STRAIGHT ADAPTER. 14 CM



CARDIOPLEGIA PRODUCTS



CSY16 STRAIGHT ADAPTER. 15 CM



CSY17 PERFUSION ADAPTER 1/4"



CSY18 PERFUSION ADAPTER 3/8"



CSY19 PERFUSION ADAPTER 1/4" WITH 3-WAY STOPCOCK



CSY20 PERFUSION ADAPTER WITH 2 FEMALE LUER LOCKS, 1 MALE LUER LOCK AND COLOR-CODED CLAMPS



CSY21 RECIRCULATION Y-ADAPTER



CSY22 MULTIPLE PERFUSION ADAPTER WITH VENT LINE



CSY23 SIX MULTIPLE PERFUSION ADAPTER



BALLOON TIP



DETAILS	
Length	30 cm
Qty/box	20 Pcs

CODE	BALLOON SIZE	TIP SIZE
OC06A	6 mm	3 mm
OC07A	7 mm	3 mm

+ Features & Benefits

- ✓ Self-inflating silicone balloon
- ✓ Female luer lock connection

BASKET TIP



OSTIAL PERFUSION CANNULAE

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

DETAILS	
Length	21 cm
Qty/box	20 Pcs

+ Features & Benefits

- ✓ Basket tip
- ✓ Different configurations
- ✓ Malleable stainless steel shaft

	TIP SIZE	CODE
WITH 45° TIP	10 Fr	OC1045
	12 Fr	OC1245
	14 Fr	OC1445
WITH 90° TIP	10 Fr	OC1090
	12 Fr	OC1290
	14 Fr	OC1490



SILICONE VENT CATHETER

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.



Features & Benefits

- ✓ With guidewire stylet
- ✓ Soft and flexible silicone tubing
- ✓ Soft tip



DETAILS	
Code	SV20
Size	20 Fr
Length	35 cm
Tip perforation	6,5 cm
Connection site	1/4"
Qty/box	20 Pcs

DETAILS	
Code	SV16
Size	16 Fr
Length	34 cm
Tip perforation	6,5 cm
Connection site	1/4"
Qty/box	20 Pcs



CLOSED TIP, RIGHT ANGLED, PREFORMED

PVC VENT CATHETER

Andocor vent catheters are medical devices used to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.

DETAILS	
Connection site	1/4"
Qty/box	20 Pcs



Features & Benefits

- ✓ Kink resistant reinforced tubing
- ✓ Easy tip orientation thanks to indexing marker

	SIZE	TIP PERFORATION	LENGTH	CODE
WITHOUT STYLET	13 Fr	3,8 cm	32 cm	PVBA13
	20 Fr	10 cm	36,5 cm	PVBA20

VENT CATHETERS



OPEN TIP, STRAIGHT



	SIZE	TIP PERFORATION	LENGTH	CODE
WITHOUT STYLET	17 Fr	2,5 cm	32 cm	PVSO170
	18 Fr	2,5 cm	32 cm	PVSO180
	17 Fr	6 cm	32 cm	PVSO171
	18 Fr	6 cm	32 cm	PVSO181
	17 Fr	9 cm	32 cm	PVSO172
	18 Fr	9 cm	32 cm	PVSO182



	SIZE	TIP PERFORATION	LENGTH	CODE
WITH PLASTIC STYLET	17 Fr	2,5 cm	32 cm	PVSOR170
	18 Fr	2,5 cm	32 cm	PVSOR180
	17 Fr	6 cm	32 cm	PVSOR171
	18 Fr	6 cm	32 cm	PVSOR181
	17 Fr	9 cm	32 cm	PVSOR172
	18 Fr	9 cm	32 cm	PVSOR182



	SIZE	TIP PERFORATION	LENGTH	CODE
WITH GUIDEWIRE STYLET	18 Fr	6 cm	32 cm	PVSOG181

CLOSED TIP, STRAIGHT



	SIZE	TIP PERFORATION	LENGTH	CODE
WITHOUT STYLET	18 Fr	2,5 cm	32 cm	PVSB180
	18 Fr	6 cm	32 cm	PVSB181
	20 Fr	6 cm	32 cm	PVSB201
	18 Fr	9 cm	32 cm	PVSB182



	SIZE	TIP PERFORATION	LENGTH	CODE
WITH GUIDEWIRE STYLET	17 Fr	2,5 cm	32 cm	PVSBG170
	18 Fr	2,5 cm	32 cm	PVSBG180
	17 Fr	6 cm	32 cm	PVSBG171
	18 Fr	6 cm	32 cm	PVSBG181
	20 Fr	6 cm	32 cm	PVSBG201
	17 Fr	9 cm	32 cm	PVSBG172
	18 Fr	9 cm	32 cm	PVSBG182



	SIZE	TIP PERFORATION	LENGTH	CODE
WITH PLASTIC STYLET	17 Fr	2,5 cm	32 cm	PVSBR170
	18 Fr	2,5 cm	32 cm	PVSBR180
	17 Fr	6 cm	32 cm	PVSBR171
	18 Fr	6 cm	32 cm	PVSBR181
	20 Fr	6 cm	32 cm	PVSBR201
	17 Fr	9 cm	32 cm	PVSBR172
	18 Fr	9 cm	32 cm	PVSBR182



RIGID SUCKER

Andocor suckers are used for the removal of excess fluid from the surgical field.

DETAILS	
Code	RS01
Length	33 cm
Qty/box	20 Pcs

+ Features & Benefits

- ✓ Metal suction tube with plastic handle
- ✓ Easy & comfortable directing



FLEXIBLE & MALLEABLE SUCKERS


Andocor suckers are used for the removal of excess fluid from the surgical field.

DETAILS	
Connection site	1/4"
Qty/box	20 Pcs

+ Features & Benefits

- ✓ Provide continuous suction
- ✓ Flexible or malleable tubing

TUBING	CODE	LENGTH	FRENCH SIZE
FLEXIBLE	FS300	36 cm	17 Fr
MALLEABLE	FS301	45 cm	17 Fr



YANKAUER SUCTION TUBE

Andocor yankauer suction tubes are intended to aspirate body waste fluid during surgery.

DETAILS	
Qty/box	20 Pcs

CONFIGURATONS	TIP	CODE
VENTED	OPEN TIP	ST04
	BULB TIP	ST02
NON-VENTED	OPEN TIP	ST03
	BULB TIP	ST01

SUCKERS



PERICARDIAL SUMPS

Andocor suckers are used for the removal of excess fluid from the surgical field.



Features & Benefits

- ✓ Smooth and flush design
- ✓ Weighted tip
- ✓ Drainage of pericardial sack

DETAILS	
Connection site	1/4"
Qty/carton	20 Pcs



PS40M | PS40ML | PS40MLH

PS40M1

METAL TIP

CODE	DESCRIPTION	TUBING LENGTH	EXTERNAL DIAMETER TIP
PS40M	METAL TIP	36 cm	6 mm
PS40ML	LONG TUBING	80 cm	6 mm
PS40MLH	LONG, HARD TUBING	80 cm	6 mm

CODE	DESCRIPTION	TUBING LENGTH	EXTERNAL DIAMETER TIP
PS40M1	LARGE METAL TIP	36 cm	7,8 mm

WIREWOUND TIP



DETAILS	
Code	PS40
Size	16 Fr
Length	36 cm
Qty/box	20 Pcs



VESSEL CANNULAE

Andocor vessel cannulae are used to help check for leaks in a harvested vein that will be used in a graft with a sterile saline solution. Additionally, the vessel cannula can be used for antegrade cardioplegia administration.

	CODE	LENGTH	QTY/BOX
TRANSPARENT, FLEXIBLE	VC3	5 cm	40 Pcs
WITH CLOSED BACK CHECK VALVE	VC3B	6,5 cm	40 Pcs
WHITE, RIGID	VC04	5 cm	40 Pcs



TUBING ORGANIZER

DETAILS	
CODE	TO10
SLOTS & TYPE	1 TUBE 1/8"
	6 TUBES 1/4"
	2 TUBES 3/8"
	1 TUBE 1/2"
QTY/BOX	10 Pcs