


Peripheral Treatment with ePTFE



FlowLine Bipore

- Bipore design with two different fibril lengths for low thrombogenicity
- An additional ePTFE wrap for enhanced suture retention and high burst strength
- Unique guideline indicates diameter and wall thickness of the graft
- Helical reinforcement for increased resistance against kinking and compression
- Excellent and pliable handling and suture behaviour
- Simple and easy removal of the spiral reinforcement

Ordering Information

		Thin Wall	Thin Wall Reinforced			Standard Wall	Standard Wall Reinforced
Length (cm)	Diameter (mm)	Catalogue No.	Catalogue No.	Length (cm)	Diameter (mm)	Catalogue No.	Catalogue No.
10	5	10TW1005N	10TW1005S	20	6	10SW2006N	–
20	6	10TW2006N	–	20	7	10SW2007N	–
20	7	10TW2007N	–	40	5	10SW4005N	–
20	8	10TW2008N	–	50	6	10SW5006N	10SW5006S
40	5	10TW4005N	10TW4005S	50	7	10SW5007N	10SW5007S
40	6	10TW4006N	–	50	8	10SW5008N	10SW5008S
50	6	10TW5006N	10TW5006S	70	5	10SW7005N	10SW7005S
50	7	10TW5007N	10TW5007S	80	6	10SW8006N	10SW8006S
50	8	10TW5008N	10TW5008S	80	7	10SW8007N	10SW8007S
70	5	10TW7005N	10TW7005S	80	8	10SW8008N	10SW8008S
80	6	10TW8006N	10TW8006S				
80	7	10TW8007N	10TW8007S				
80	8	10TW8008N	10TW8008S				

ARTIVION™

Learn more at [artivion.com](https://www.artivion.com)

All products and indications are not available/approved in all markets. All trademarks are owned by Artivion, Inc. or its subsidiaries. On-X Life Technologies, Inc. Jotec GmbH and Ascyrus Medical GmbH are wholly owned subsidiaries of Artivion, Inc. © 2022 Artivion, Inc. All rights reserved.

Artivion, Inc.

1655 Roberts Blvd., NW
Kennesaw, GA 30144 USA
Phone: +1-888-427-9654
Fax: +770-590-3753

E-mail: inquiries@artivion.com

For contact information by region, please visit www.artivion.com/contact

 JOTEC GmbH, Lotzenäcker 23, 72379 Hechingen, Germany

JT-BR-0100200-EN V03 12/2022