

Hypogastric Artery Matters



ARTIVION™

E-liac®
Stent Graft System

Ordering Information

E-lic Stent Graft System

01 Aorto-Iliac Aneurysms



Catalog Number	Proximal Ø (mm)	Distal Ø (mm)	Total Length (mm)	Proximal Length (mm)	Distal Length (mm)	OD Delivery System (F/mm)
72IB1410L53L44	14	10	97	53	44	18/6.25
72IB1412L53L44	14	12	97	53	44	18/6.25
72IB1414L53L44	14	14	97	53	44	18/6.25

Special sizes on demand only

72IB1410L53L56	14	10	109	53	56	18/6.25
72IB1412L53L56	14	12	109	53	56	18/6.25
72IB1414L53L56	14	14	109	53	56	18/6.25

02 Isolated Iliac Aneurysms



Catalog Number	Proximal Ø (mm)	Distal Ø (mm)	Total Length (mm)	Proximal Length (mm)	Distal Length (mm)	OD Delivery System (F/mm)
72IB1410L65L44	14	10	109	65	44	18/6.25
72IB1412L65L44	14	12	109	65	44	18/6.25
72IB1414L65L44	14	14	109	65	44	18/6.25
72IB1612L65L56	16	12	121	65	56	18/6.25
72IB1812L65L56	18	12	121	65	56	18/6.25

Special sizes on demand only

72IB1412L65L56	14	12	121	65	56	18/6.25
72IB1414L65L56	14	14	121	65	56	18/6.25
72IB1610L65L56	16	10	121	65	56	18/6.25
72IB1614L65L56	16	14	121	65	56	18/6.25
72IB1814L65L56	18	14	121	65	56	18/6.25

ARTIVION™

Learn more at artivion.com

All products and indications are not available/approved in all markets. All trademarks are owned by Artivion, Inc. or its subsidiaries. On-X Life Technologies, Inc. Jotec GmbH and Ascyrus Medical GmbH are wholly owned subsidiaries of Artivion, Inc. Artivion and E-lic are Trademarks of Artivion, Inc. or its subsidiaries. E-lic is registered in Germany and/or other countries. © 2022 Artivion, Inc. All rights reserved.

Artivion, Inc.

1655 Roberts Blvd., NW
Kennesaw, GA 30144 USA
Phone: +1-888-427-9654
Fax: +770-590-3753

E-mail: inquiries@artivion.com

For contact information by region, please visit www.artivion.com/contact

 JOTEC GmbH, Lotzenäcker 23, 72379 Hechingen, Germany

JT-RS-0720200-EN V01 01/2023