

## Disposable Incision Protection Sleeve ( Disposable Wound Retractors )

### A TYPE



Spec.	Outer hoop Size(mm)	Inner hoop Size(mm)	Channel size(Outer)(mm)	Channel Height(mm)
LMQK-60/70-60/150(B)	60	70	60	Up to 150
LMQK-80/90-80/150(B)	80	90	80	Up to 150
LMQK-120/130-120/250(B)	120	130	120	Up to 250
LMQK-150/160-150/250(B)	150	160	150	Up to 250
LMQK-180/190-180/250(B)	180	190	180	Up to 250
LMQK-200/210-200/250(B)	200	210	200	Up to 250
LMQK-220/230-220/250(B)	220	230	220	Up to 250
LMQK-250/260-250/250(B)	250	260	250	Up to 250
LMQK-270/280-270/250(B)	270	280	270	Up to 250
LMQK-320/330-320/250(B)	320	330	320	Up to 250

### B TYPE



Spec.	Outer hoop Size(mm)	Inner hoop Size(mm)	Channel size(Outer)(mm)	Channel Height(mm)
LMQK-40/30-15/25 (D)	40	30	15	25
LMQK-40/30-15/30 (D)	40	30	15	30
LMQK-40/30-15/35 (D)	40	30	15	35
LMQK-50/40-25/25 (D)	50	40	25	25
LMQK-50/40-25/30 (D)	50	40	25	30
LMQK-50/40-25/35 (D)	50	40	25	35
LMQK-70/70-35/25 (D)	70	70	35	25
LMQK-70/70-35/30 (D)	70	70	35	30
LMQK-70/70-35/35 (D)	70	70	35	35
LMQK-100/80-60/25 (D)	100	80	60	25
LMQK-100/80-60/40 (D)	100	80	60	40
LMQK-120/100-60/25 (D)	120	100	60	25
LMQK-120/100-60/40 (D)	120	100	60	40

### C TYPE



Model	Outer ring	Inner ring
LMQK50/120	120	50

#### Features include:

- ※ LMQK-A Made of PE and TPU on standard GB / T11115
- ※ LMQK-A Suit for different size operation in endoscope and small incision operation.
- ※ LMQK-B Made of silica gel and 12Cr18Ni9 SUS wire.
- ※ LMQK-B Suit for thoracoscope, laparoscope, breast, thyroid gland and microinvasive operation.